



Reconditioned Equipment

Many customers are looking to replace equipment but have limited budgets and for this reason Prescott Instruments have become well established to rubber companies in over 20 countries.

In our own modern workshop, instruments are stripped down completely before being resprayed and rebuilt with new components including modern digital PID temperature control systems, new dies & rotors and new pneumatic systems. In addition many instruments can be supplied with a new PC and software system giving the customer a piece of equipment that will give good service for many years to come at a fraction of the cost of a new unit.

All Prescott Instruments refurbished equipment carries a full warranty.

Enjoy the benefits of quality laboratory equipment at a fraction of the cost of new.

In addition to supplying reconditioned/upgraded equipment from stock, we can also undertake to refurbish customers existing equipment either on site or at our premises, dependant upon specific requirements.

If customers do not want the inconvenience of being without their testing equipment during refurbishment, part exchange is an option which may be considered. Loan instruments are also available.



The variety of equipment we refurbish includes:

- Monsanto R100 R100s Rheometer
- Monsanto 1500s Rheometer
- Negretti Mooney Mk3, Mk8
- Alpha ODR 2000
- Alpha MDR 2000
- Alpha MV 2000

Please note reconditioned equipment is subject to availability.

Notice:

Although the information and recommendations set forth herein (hereinafter "Information") are presented in good faith and believed to be correct as of the date hereof. Prescott Instruments Ltd. make no representation or warranties as to the completeness or accuracy thereof. Information is supplied upon the conditions that the persons receiving same will make their own determination as to the suitability for their purposes prior to use. In no event will Prescott Instruments Ltd. be responsible for damages of any nature whatsoever resulting from the use or reliance upon information for the product, equipment or system to which information refers. Nothing contained herein is to be construed as recommendation to use any product, equipment, system, process or formulation in conflict with any patent, and Prescott Instruments Ltd. makes no representation or warranty, express or implied, that the use thereof will not infringe any patent, except for the limited warranty set forth in Prescott Instruments Ltd. standard sales contracts for its equipment and services.

PRESCOTT INSTRUMENTS LTD. MAKES NO REPRESENTATIONS OR WARRANTIES WHETHER STATUTORY, EXPRESS OR IMPLIED, OF MERCHANTABILITY FITNESS FOR A PARTICULAR PURPOSE, OR OF ANY OTHER NATURE WITH RESPECT TO THE INFORMATION OR PRODUCT, EQUIPMENT OR SYSTEM TO WHICH INFORMATION REFERS.



**PRESCOTT
INSTRUMENTS LTD**

Prescott Instruments Ltd. Unit F, Northway Trading Estate,
Northway Lane, Tewkesbury, Gloucester GL20 8JH U.K.

Tel: +44 (0) 1684 274300 Fax: +44 (0) 1684 293223

Email: enquiries@prescott-instruments.com

Website: www.prescott-instruments.com

