


Increasingly, laboratories are looking to integrate their test information into a more complete data network combining test data with many other measurements taken throughout the production process and Lab Line is ideally suited to this application.

Prescott Instruments' Lab Line system is a fully integrated approach, offering a network system connecting any type of rubber testing instrument to a centralised database.

The database type can be chosen to suit the size of the installation, Microsoft Access for up to 4 instruments and SQL or IBM DB2 can be used for larger installations.

Labline software can be used to connect a number of separate laboratories together and streamline the flow of data enabling reports to be generated.



Integration with other technologies such as SAP is also a feature of this very flexible package.

Real time displays with auto scaling of the graph axis gives the best possible presentation of testing information of test data combined with automatic pass / fail at the end of each test with fully programmable limits.

The system can accommodate any number of test instruments with all of the data held in a secure Access database.

